

Por infrarrojos | **Infrared thermometer**

Marca|**Trademark:** Nahita

El termómetro clínico digital por infrarrojos de Nahita permite medir la temperatura corporal en el oído o la frente de forma fiable y en cuestión de segundos gracias al sensor localizado en la punta del termómetro. Gracias a su diseño, el manejo del termómetro es sencillo mediante un único botón y no es necesario el uso de protectores desechables del sensor.

Está equipado con un capuchón para la determinación de la temperatura corporal en la frente y un práctico soporte para tener el termómetro siempre a mano para su utilización.

Con memoria de la última lectura realizada que se visualiza cada vez que se enciende el termómetro y alarma sonora en caso de fiebre (temperatura corporal por encima de 38 °C).

Puede ser también programado para mostrar la hora y temperatura ambiental.

Digital clinical infrared thermometer NAHITA measures body temperature in the ear or forehead reliably and in a matter of seconds thanks to the sensor situated at the tip of the thermometer. Thanks to its design, the thermometer use is easy by a unique button and the use of disposable sensor covers is not necessary. It is equipped with a probe cap for forehead temperature determination and a useful support to have always the thermometer handy.

With memory of the last read performed that is displayed every time the thermometer is turned on and an acoustic alarm in case of fever (corporal temperature higher than 38 °C). It can also be programmed to show the hour and environmental temperature



Referencia / Code	72752500
Rango temperatura corporal / Corporal temperature range	32-42.9 °C / 89.6-109 °F
Exactitud / Accuracy	± 0.2 °C (35.5-41.0 °C) ± 0.3 °C (32.0-35.4 °C ó 41.1-42.9 °C)
Rango temperatura ambiental / Environmental temperature range	10-40 °C
Precisión / Precision	0.1 °C / 0.1 °F
Apagado automático / Auto off function	Sí / Yes
Batería / Battery	1 x 3 V, litio, tipo CR2032 / 1 x 3 V, lithium, CR2032 type